


Product Datasheet

Rat ATG5 (autophagy protein 5) ELISA Kit (orb1100797)

Description	Rat ATG5 (autophagy protein 5) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	ATG5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q3MQ06
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat ATG5
(autophagy
protein 5)
ELISA Kit