


Product Datasheet

Rat IL15RA (Interleukin-15 receptor subunit alpha) ELISA Kit (orb1100800)

Description	Rat IL15RA (Interleukin-15 receptor subunit alpha) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	IL15RA
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	F1M293
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat IL15RA
(Interleukin-
15 receptor
subu...