


Product Datasheet

Mouse Ralbp1 (RalA-binding protein 1) ELISA Kit (orb1100809)

Description	Mouse Ralbp1 (RalA-binding protein 1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Ralbp1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q62172
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Ralbp1
(RalA-
binding
protein 1)
EL...