


Product Datasheet

Mouse Pglyrp1 (Peptidoglycan recognition protein 1) ELISA Kit (orb1100810)

Description	Mouse Pglyrp1 (Peptidoglycan recognition protein 1) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Pglyrp1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O88593
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse Pglyrp1
(Peptidoglycan
recognition...