


Product Datasheet

Mouse Hmgcs2 (Hydroxymethylglutaryl-CoA synthase, mitochondrial) ELISA Kit (orb1100811)

Description	Mouse Hmgcs2 (Hydroxymethylglutaryl-CoA synthase, mitochondrial) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Hmgcs2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P54869
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Hmgcs2
(Hydroxymethylglutaryl-CoA ...