


Product Datasheet

Mouse Hn1l (Hematological and neurological expressed 1-like protein) ELISA Kit (orb1100813)

Description	Mouse Hn1l (Hematological and neurological expressed 1-like protein) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Hn1l
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q6PGH2
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse Hn1l
(Hematological
and
neurologic...