


Product Datasheet

Mouse Enpep (Glutamyl aminopeptidase) ELISA Kit (orb1100814)

Description	Mouse Enpep (Glutamyl aminopeptidase) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Enpep
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P16406
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse Enpep
(Glutamyl
aminopeptidase)
EL...