


Product Datasheet

Mouse Ctps1 (CTP synthase 1) ELISA Kit (orb1100818)

Description	Mouse Ctps1 (CTP synthase 1) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Ctps1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P70698
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Ctps1
(CTP
synthase
1) ELISA
Kit