


## Product Datasheet

### Mouse C3AR1 (C3a anaphylatoxin chemotactic receptor) ELISA Kit (orb1100819)

|                        |                                                                |
|------------------------|----------------------------------------------------------------|
| <b>Description</b>     | Mouse C3AR1 (C3a anaphylatoxin chemotactic receptor) ELISA Kit |
| <b>Reactivity</b>      | Mouse                                                          |
| <b>Range</b>           | 78.125-5000pg/ml                                               |
| <b>Target</b>          | C3AR1                                                          |
| <b>Storage</b>         | 2-8 °C for 6 months                                            |
| <b>Note</b>            | For research use only                                          |
| <b>Uniprot ID</b>      | <a href="#">O09047</a>                                         |
| <b>Sensitivity</b>     | 46.875pg/ml                                                    |
| <b>Expiration Date</b> | Please enquire.                                                |



Mouse C3AR1  
(C3a  
anaphylatoxin  
chemotact...