


Product Datasheet

Mouse Cerk (Ceramide kinase) ELISA KIT (orb1100820)

Description	Mouse Cerk (Ceramide kinase) ELISA KIT
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Cerk
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8K4Q7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Cerk
(Ceramide
kinase)
ELISA KIT