


Product Datasheet

Mouse Adora1 (Adenosine receptor A1) ELISA Kit (orb1100825)

Description	Mouse Adora1 (Adenosine receptor A1) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Adora1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q60612
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Adora1
(Adenosine
receptor
A1) ELI...