


Product Datasheet

Mouse MMP14 (Matrix Metalloproteinase-14) ELISA kit (orb1100829)

Description	Mouse MMP14 (Matrix Metalloproteinase-14) ELISA kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	MMP14
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P53690
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse MMP14
(Matrix
Metalloproteinase-
14...