


Product Datasheet

Mouse ANGPT4(Angiopoietin-4) ELISA Kit (orb1100831)

Description	Mouse ANGPT4(Angiopoietin-4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	ANGPT4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9WVH6
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
ANGPT4(Angiopoietin-
4) ELISA Kit