


Product Datasheet

Mouse Nr3c2(Mineralocorticoid receptor) ELISA Kit (orb1100833)

Description	Mouse Nr3c2(Mineralocorticoid receptor) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Nr3c2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8VII8
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Nr3c2(Mineralocorticoid
receptor) ...