


Product Datasheet

Mouse Phex (Phosphate-regulating neutral endopeptidase PHEX)ELISA Kit (orb1100844)

Description	Mouse Phex (Phosphate-regulating neutral endopeptidase PHEX)ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Phex
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P70669
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse Phex
(Phosphate-
regulating
neutral...