


Product Datasheet

Mouse Npc1 (NPC intracellular cholesterol transporter 1) ELISA Kit (orb1100845)

Description	Mouse Npc1 (NPC intracellular cholesterol transporter 1) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Npc1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O35604
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse Npc1
(NPC
intracellular
cholestero...