


Product Datasheet

Mouse LAMP2(Lysosome-associated membrane glycoprotein 2) ELISA Kit (orb1100846)

Description	Mouse LAMP2(Lysosome-associated membrane glycoprotein 2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	LAMP2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P17047
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
LAMP2(Lysosome-associated
membrane...