


Product Datasheet

Mouse CHOP(C/EBP Homologous Protein) Elisa Kit (orb1100847)

Description	Mouse CHOP(C/EBP Homologous Protein) Elisa Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	CHOP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P35639
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
CHOP(C/EBP
Homologous
Protein)
Eli...