


Product Datasheet

Mouse Ccl9 / MIP-1 Gamma (C-C motif chemokine 9) ELISA Kit (orb1100849)

Description	Mouse Ccl9 / MIP-1 Gamma (C-C motif chemokine 9) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Ccl9
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P51670
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Ccl9 /
MIP-1
Gamma
(C-C motif
chem...