


Product Datasheet

Mouse AQP9(Aquaporin-9) ELISA Kit (orb1100850)

Description	Mouse AQP9(Aquaporin-9) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Aqp9
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9JJJ3
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
AQP9(Aquaporin-
9) ELISA Kit