


Product Datasheet

Mouse Slc22a12 (Solute carrier family 22 member 12) ELISA Kit (orb1100851)

Description	Mouse Slc22a12 (Solute carrier family 22 member 12) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Slc22a12
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8CFZ5
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Slc22a12
(Solute
carrier
family
22...