


Product Datasheet

Mouse IL-19 (Interleukin-19) ELISA Kit (orb1100863)

Description	Mouse IL-19 (Interleukin-19) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	IL-19
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8CJ70
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse IL-19
(Interleukin-
19) ELISA
Kit