


Product Datasheet

Mouse SFRP5 (Secreted frizzled-related protein 5) ELISA Kit (orb1100866)

Description	Mouse SFRP5 (Secreted frizzled-related protein 5) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	SFRP5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9WU66
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
SFRP5
(Secreted
frizzled-
related
p...