


Product Datasheet

Mouse D-LDH(D-Lactate Dehydrogenase) ELISA Kit (orb1100873)

Description	Mouse D-LDH(D-Lactate Dehydrogenase) ELISA Kit
Reactivity	Mouse
Range	4.688-300mIU/ml
Target	D-LDH
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7TNG8
Sensitivity	2.813mIU/ml
Expiration Date	Please enquire.



Mouse D-LDH(D-Lactate Dehydrogenase) ELI...