


Product Datasheet

Hamster M-CSF (Macrophage Colony-Stimulating Factor)ELISA Kit (orb1100882)

Description	Hamster M-CSF (Macrophage Colony-Stimulating Factor)ELISA Kit
Reactivity	Hamster
Range	15.625-1000pg/ml
Target	M-CSF
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A0A1U7QLU2
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Hamster M-
CSF
(Macrophage
Colony-
Stimula...