


Product Datasheet

Human MAGEA1 (Melanoma-associated antigen 1) ELISA Kit (orb1100885)

Description	Human MAGEA1 (Melanoma-associated antigen 1) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	MAGEA1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P43355
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
MAGEA1
(Melanoma-
associated
antige...