


Product Datasheet

Human LSP1(Lymphocyte-specific protein 1) ELISA Kit (orb1100886)

Description	Human LSP1(Lymphocyte-specific protein 1) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	LSP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P33241
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
LSP1(Lymphocyte-specific protein
1...