


Product Datasheet

Human LBR(Delta(14)-sterol reductase LBR) ELISA Kit (orb1100887)

Description	Human LBR(Delta(14)-sterol reductase LBR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	LBR
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q14739
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
LBR(Delta(14)-
sterol
reductase
LBR...