


Product Datasheet

Human IRF4 (Interferon Regulatory Factor 4)ELISA Kit (orb1100891)

Description	Human IRF4 (Interferon Regulatory Factor 4)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	IRF4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15306
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
IRF4
(Interferon
Regulatory
Factor...