


Product Datasheet

Human INO80 (Chromatin-remodeling ATPase INO80) ELISA Kit (orb1100895)

Description	Human INO80 (Chromatin-remodeling ATPase INO80) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	INO80
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9ULG1
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
INO80
(Chromatin-remodeling
ATPase...)