


Product Datasheet

Human HSD3B2 (3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 2) ELISA Kit (orb1100898)

Description	Human HSD3B2 (3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 2) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	HSD3B2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P26439
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
HSD3B2 (3
beta-
hydroxysteroid
dehy...