


Product Datasheet

Human HOMER1(Homer protein homolog 1) ELISA Kit (orb1100899)

Description	Human HOMER1(Homer protein homolog 1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	HOMER1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q86YM7
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
HOMER1(Homer
protein
homolog 1)
EL...