


Product Datasheet

Human HIVEP1(Zinc finger protein 40) ELISA Kit (orb1100901)

Description	Human HIVEP1(Zinc finger protein 40) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	HIVEP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P15822
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
HIVEP1(Zinc
finger
protein 40)
ELI...