


Product Datasheet

Human GULP1 (PTB domain-containing engulfment adapter protein 1) ELISA Kit (orb1100902)

Description	Human GULP1 (PTB domain-containing engulfment adapter protein 1) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	GULP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UBP9
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
GULP1
(PTB
domain-
containing
engul...