


Product Datasheet

Human GPX4 (Phospholipid hydroperoxide glutathione peroxidase)ELISA Kit (orb1100904)

Description	Human GPX4 (Phospholipid hydroperoxide glutathione peroxidase)ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	GPX4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P36969
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human GPX4
(Phospholipid
hydroperoxide
g...