


## Product Datasheet

### Human GAGE7(G antigen 7) ELISA Kit (orb1100909)

|                        |                                    |
|------------------------|------------------------------------|
| <b>Description</b>     | Human GAGE7(G antigen 7) ELISA Kit |
| <b>Reactivity</b>      | Human                              |
| <b>Range</b>           | 0.156-10ng/ml                      |
| <b>Target</b>          | GAGE7                              |
| <b>Storage</b>         | 2-8 °C for 6 months                |
| <b>Note</b>            | For research use only              |
| <b>Uniprot ID</b>      | <b>076087</b>                      |
| <b>Sensitivity</b>     | 0.094ng/ml                         |
| <b>Expiration Date</b> | Please enquire.                    |



Human  
GAGE7(G  
antigen  
7) ELISA  
Kit