


Product Datasheet

Human GABARAPL1 (Gamma-aminobutyric acid receptor-associated protein-like 1) ELISA Kit (orb1100910)

Description	Human GABARAPL1 (Gamma-aminobutyric acid receptor-associated protein-like 1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	GABARAPL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H0R8
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
GABARAPL1
(Gamma-aminobutyric acid...)