


Product Datasheet

Human GABARAP(Gamma-aminobutyric acid receptor-associated protein) ELISA Kit (orb1100911)

Description	Human GABARAP(Gamma-aminobutyric acid receptor-associated protein) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	GABARAP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O95166
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
GABARAP(Gamma-aminobutyric acid
receptor-associated protein) ELISA Kit