


Product Datasheet

Human Cyclin A1(Cyclin-A1) ELISA Kit (orb1100919)

Description	Human Cyclin A1(Cyclin-A1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	Cyclin A1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P78396
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
Cyclin
A1(Cyclin-
A1) ELISA
Kit