


Product Datasheet

Human CXXC5(CXXC-Type Zinc Finger Protein 5) ELISA Kit (orb1100920)

Description	Human CXXC5(CXXC-Type Zinc Finger Protein 5) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CXXC5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7LFL8
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CXXC5(CXXC-
Type Zinc
Finger
Protei...