


Product Datasheet

Human CREB3L2(Cyclic AMP-responsive element-binding protein 3-like protein 2) ELISA Kit (orb1100922)

Description	Human CREB3L2(Cyclic AMP-responsive element-binding protein 3-like protein 2) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	CREB3L2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q70SY1
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
CREB3L2(Cyclic
AMP-
responsive
elem...