


Product Datasheet

Human CAMK1D(Calcium/calmodulin-dependent protein kinase type 1D) ELISA Kit (orb1100925)

Description	Human CAMK1D(Calcium/calmodulin-dependent protein kinase type 1D) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	CAMK1D
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8IU85
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
CAMK1D(Calcium/calmodulin-
dependen...