


## Product Datasheet

### Human ARPC5(Actin-related protein 2/3 complex subunit 5) ELISA Kit (orb1100933)

|                        |                                                                    |
|------------------------|--------------------------------------------------------------------|
| <b>Description</b>     | Human ARPC5(Actin-related protein 2/3 complex subunit 5) ELISA Kit |
| <b>Reactivity</b>      | Human                                                              |
| <b>Range</b>           | 78.125-5000pg/ml                                                   |
| <b>Target</b>          | ARPC5                                                              |
| <b>Storage</b>         | 2-8 °C for 6 months                                                |
| <b>Note</b>            | For research use only                                              |
| <b>Uniprot ID</b>      | <b>O15511</b>                                                      |
| <b>Sensitivity</b>     | 46.875pg/ml                                                        |
| <b>Expiration Date</b> | Please enquire.                                                    |



Human  
ARPC5(Actin-  
related  
protein 2/3  
co...