


Product Datasheet

Human ADH6(Alcohol dehydrogenase 6) ELISA Kit (orb1100934)

Description	Human ADH6(Alcohol dehydrogenase 6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	ADH6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P28332
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
ADH6(Alcohol
dehydrogenase
6) ELIS...