


Product Datasheet

Human Anti-2019 nCoV(S) IgA ELISA Kit (orb1100940)

Description	Human Anti-2019 nCoV(S) IgA ELISA Kit
Reactivity	Human
Range	0.391-25ng/ml
Target	Anti-2019 nCoV(S) IgA
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.234ng/ml
Expiration Date	Please enquire.



Human
Anti-2019
nCoV(S)
IgA ELISA
Kit