


Product Datasheet

Human anti-2019 nCoV(N+S) IgM ELISA Kit (orb1100942)

Description	Human anti-2019 nCoV(N+S) IgM ELISA Kit
Reactivity	Human
Range	0.781-50ng/ml
Target	anti-2019 nCoV(N+S) IgM
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Indirect Elisa
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Human
anti-2019
nCoV(N+S)
IgM ELISA
Kit