


Product Datasheet

Human GRM5 (Metabotropic glutamate receptor 5)ELISA Kit (orb1100945)

Description	Human GRM5 (Metabotropic glutamate receptor 5)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	GRM5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P41594
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human GRM5
(Metabotropic
glutamate
recep...