


Product Datasheet

Human GABRB2 (Gamma-aminobutyric acid receptor subunit beta-2)ELISA Kit (orb1100946)

Description	Human GABRB2 (Gamma-aminobutyric acid receptor subunit beta-2)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	GABRB2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P47870
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
GABRB2
(Gamma-aminobutyric acid re...