


Product Datasheet

Human M2BPGI(Mac-2 Binding Protein Glycosylation Isomer) ELISA Kit (orb1100949)

Description	Human M2BPGI(Mac-2 Binding Protein Glycosylation Isomer) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	M2BPGI
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
M2BPGI(Mac-
2 Binding
Protein
Glyco...