


Product Datasheet

Human SHANK3(SH3 and multiple ankyrin repeat domains protein 3)ELISA Kit (orb1100964)

Description	Human SHANK3(SH3 and multiple ankyrin repeat domains protein 3)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SHANK3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9BYB0
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human SHANK3(SH3 and multiple ankyrin re...